



TECHNOLOGY CENTER

GC Analysis and Techniques



CTE GLOBAL, INC.

Leading to a brighter tomorrow



Presentation Outline

1) Basic overview of GC

- Separation Principle
- Carrier gas
- Maintenance and consumables


2) ASTM D5501

- Running parameters
- Calibrations
- Calculations


3) CTE Process Monitoring Services

- Profiling fusel production
- Amino acid profile

GC Basics – Separation Principle



Separation Of Compounds



- Stationary phase and mobile phase are fixed throughout analysis
- Carrier gas flow rate and oven temperature are the only variables changing throughout analysis
- The separation of compounds is heavily dependent on its boiling point and its interaction with the stationary phase
- Compounds with a lower boiling point will travel through the column faster and elute first
- If boiling points are similar, then the compounds will elute with respect to their polarity – weaker interaction with the stationary phase will elute first

Separation Of Compounds

- Separation of compounds in a sample matrix depends mainly on two factors:
 - Column temperature, each compound elutes at different times based on its boiling point
 - The polarity of a components in a sample vs polarity of stationary phase

- But also depends on:
 - Carrier gas
 - Column length
 - Injection volume
 - + more

This will cause the compounds to elute at different times (retention time) and flow to the detector

Carrier Gas

To achieve good separation and reproducible chromatograms:

1. Carrier gas supply needs to be high-purity
 - Min. purity 99.95% - Reduces base line noise
2. Carrier gas flow needs to be constant
 - Using pressure and flow regulators
3. Carrier gas needs to be inert
 - so, it won't interfere with the separation and detection

Helium, hydrogen and nitrogen are some of the commonly-used carrier gases.

Carrier Gas

Carrier Gas	Advantages	Disadvantages
Helium	<ul style="list-style-type: none">• Inert (safe) and non-flammable• Gives high resolution	<ul style="list-style-type: none">• Expensive, not easily available
Hydrogen	<ul style="list-style-type: none">• High diffusivity and linear velocities• Gets good separation efficiencies• Short analysis and run time (results in cheap operational cost)	<ul style="list-style-type: none">• Flammable• Not completely inert (e.g. reacts with some compounds at high temperature)
Nitrogen	<ul style="list-style-type: none">• Cheap and easily available	<ul style="list-style-type: none">• Not suited for use in temperature-programmed GC analysis• Lower or poor separation resolution• Long analysis and run time


Carrier Gas – Alternatives to Helium

If using an alternative, you will need to consider:

- Safety (hydrogen requires additional safety measures)
- Purity
- Any hardware changes
- Parameter conditions

Carrier Gas

Optimize helium conservation methods through Gas Saver mode and Ecology Mode on Shimadzu GC



Equipment Maintenance/Consumables


- Should have yearly scheduled PM's
- Routinely maintain consumables

1. Gas
 - Pure, make sure you have inventory
2. Syringe
 - Both for manual and autosample.
3. Septa (after 100 injections)
 - Low bleed and optimum sealing



4. Glass Liner (after 500 injections)

- Elect correct one for improved performance, vaporization



5. O-rings

- Low bleed and optimum sealing



6. Ferrules

- Ensures best connection is achieved



7. Columns


- Ensure correct installation, conditioning



GC Peak Troubleshooting

Column Bleed

- Conditioning
- Contamination
- Leak




Reduced Size

- Clogged syringe
- Leak
- Split ratio
- Temp




Tailing

- Contamination
- Column Installation/cut
- Temp
- Split ratio
- Inlet overloading




Noise

- Column installation
- Leak
- Contamination




Flap Top

- Detector overloaded




Fronting

- Column overloaded



Peak Overlap

- Column dimensions
- Column trimming
- Damaged column stationary phase



ASTM D5501

- At ethanol plants, GC analysis is used to determine the amount of ethanol, methanol and denaturant is found in fuel blends
- The standard GC test method used is ASTM D5501
Determination of Ethanol and Methanol Content in Fuels Containing Greater than 20 % Ethanol by Gas Chromatography
- Ensuring the GC is running at optimal conditions for this method is crucial for reproducibility and accuracy

ASTM D5501 Running Parameters

- Need a GC capable of operating at the conditions listed in the adjacent table
- Any column with good resolution and selectivity can be used
- Ensure consumables are within usage limits
- Accurate sample injection and split ratios is crucial to the precision and accuracy
- Carrier gas + detector gases should have a min purity of 99.95%


TABLE 1 Typical Operating Conditions

Column Temperature Program		
Column length	100 m	150 m
Initial temperature	15 °C	60 °C
Initial hold time	12 min	15 min
Program rate	30 °C/min	30 °C/min
Final temperature	250 °C	250 °C
Final hold time	19 min	23 min
Injector		
Temperature	300 °C	
Split ratio	200:1	
Sample size	0.1 µL to 0.5 µL	
Detector		
Type	Flame ionization	
Temperature	300 °C	
Fuel gas	Hydrogen (30 mL/min)	
Oxidizing gas	Air (300 mL/min)	
Make-up gas	Helium or Nitrogen (30 mL/min)	
Date rate	20 Hz	
Carrier Gas		
Type	Helium or Hydrogen ^A	
Average linear velocity	21 cm/s to 24 cm/s (constant flow)	

^A Use of hydrogen carrier gas requires additional safety considerations.


ASTM D5501 - Qualitative Analysis

- Inject ethanol and method standards
- The retention time of each standard under set parameter conditions will be used to identify the eluting compounds from your sample
- When analyzed under the same conditions, the same compound elutes at the same time



ASTM D5501 – Calibration Linearity

- Run calibration standards that cover the expected range of ethanol and methanol
- Analyze the peak areas vs concentrations and ensure there is a linear regression with a minimum $r^2=0.995$
- When peak area = 0, ethanol concentration is +/-3% by mass
- If the conditions above are not met, troubleshoot until they are
 - Increase split ratio
 - Reduce injection volume
 - Correct peak integration of ethanol and methanol




ASTM D5501 Quantitative Analysis

- Absolute Calibration
 - Single point
 - Calculate the mass response factor (MRF) of the compounds in each calibration standard
 - The average MRF is used for the calibration
 - Multi-point
 - Generate individual calibration curves for each compound in the GC software
- Relative Calibration
 - Calculate the MRF of ethanol, methanol and heptane
 - Then, calculate the mass relative response factor (MRRF) of methanol and ethanol relative to heptane
 - The average MRRF is used for the calibration

ASTM D5501 Quality Control


- Run a quality control (QC) sample daily along side samples
- Compile the data and build a control chart
- Monitor the control chart for any outliers
- Ensures precision and accuracy




ASTM D5501 Calculations

- From the standard chromatograms, obtain the peak area and calculate the response factor using mass %
- Calculate the MRRF
- From the denatured ethanol chromatographs, obtain the peak area
- Calculate for corrected peak area then normalize mass % and correct for moisture (ASTM D1364)
- Report results in volume % using compound density and sample density (ASTM D1298 or D4052)



Process Monitoring Services



GC-MS



Mass Spec Results




Type	m/z	Area	Set%	Act. %	Ref. Band
Target	43.00	6035490	100.00	100.00	---
Ref. Ion1	41.00	4642068	81.64	76.91	30
Ref. Ion2	31.00	3390469	79.87	56.18	30
Ref. Ion3	33.00	3362711	67.20	55.72	30

- Ions are separated based on their mass-to-charge ratio (m/z)
- The relative intensity of each ion is measured and then recorded to produce a mass spec
- The result displays the relative ion intensity against their m/z
- Compounds are identified by a fingerprint of mass fragments measured by mass spec

NIST Library

- Each fragment is logged into an online database (NIST library) that can identify the compound
 - The library will produce a % match, the higher the percentage, the more closely it relates to the compound being identified
- In this example, the mass fragment fingerprint matches up very well with the expected fingerprint for isobutanol





CTE TECHNOLOGY CENTER

Profiling Fusel Production




CTE GLOBAL, INC.


Leading to a brighter tomorrow

The Ehrlich Pathway

- Fusels are normal metabolites produced by yeast and bacteria, but can be about 15 times more inhibitory to yeast than ethanol
- Depending on several factors, inhibition may be evident even in the 100-500 ppm range
- The carbon skeletons of some amino acids are not incorporated into central metabolism, and are instead metabolized to fusels through the Ehrlich Pathway
- The Ehrlich pathway describes the harvesting of the nitrogen through transamination, followed by decarboxylation and oxidation to produce a fusel alcohol



Fusel Chromatograph




Compounds Identified:

- Acetaldehyde
- Methanol
- Isobutanol
- Butanol
- 3-methyl-1-butanol (Isoamyl alcohol)
- 2-methyl-1-butanol
- Pentanol
- 2,3-Butanediol
- Phenylethyl Alcohol
- Acetal
 - Can add any additional compounds upon request.

Excellent linearity

Isobutanol




$Y = 1.056087X + 9.391085e-002$
 $R^2 = 0.999737$
 $R = 0.9998685$

Internal Standard
Curve: Linear
Origin: Not Forced
Weighting Method: None

Mean RF : 1.105766
RF SD : 0.0262968
RF %RSD : 2.378150

Phenylethyl Alcohol




$Y = 2.803818X - 1.691603e-002$
 $R^2 = 0.9997972$
 $R = 0.9998986$

Internal Standard
Curve: Linear
Origin: Not Forced
Weighting Method: None

Mean RF : 2.715642
RF SD : 0.0464885
RF %RSD : 1.711878

Fusels & Organic Impurities Profile in Fermentation



Summary: Fusel Analysis

- Assess fusels production throughout fermentation
- Optimize fusels removal during distillation
- Evaluate the potential for fusels recycling in your process
- Ultimately, prevent fusels inhibition at your plant



CTE TECHNOLOGY CENTER

Amino Acid Profile For Best Performance




CTE GLOBAL, INC.


Leading to a brighter tomorrow

Nitrogen metabolism in yeast

- Research literature shows that not all nitrogen sources are the same
- Each nitrogen source can be categorized based on their ability to support growth and activate different metabolic pathways
- The diagram on the right shows the preferred nitrogen sources in green, intermediately preferred in black, and non-preferred amino acids in red



Amino acid analysis by GC-MS



- GC-MS analysis of derivatized samples resolves 21 primary amino acids and many others
- Sample prep in 8 minutes and GC-MS runs are 30 minutes
- Understand the kinetics of amino acid production and consumption in fermentations/liquefactions treated with protease
- Important insights of effects of protease dosing and urea/ammonia reductions on fusel alcohol production

Amino acid profile: Optimize performance

*Preferred amino acids: alanine, asparagine, aspartate, glutamine, glutamate, and serine

*Intermediately preferred amino acids: glycine, ornithine, phenylalanine, and valine





*Non-preferred amino acids: methionine, isoleucine, leucine, threonine, tryptophan, and tyrosine

*Not involved in growth support: histidine and lysine

Amino Acid	Table 1. Amino acid profile at each sample point (nmol/mL)										
	Liq	1hr	5.5hr	11.5hr	17.5hr	23 hr	29hr	35hr	41hr	46.5hr	TS
ALA*	1189	1539	1610	29	52	190	307	543	755	1255	2375
ASN*	1378	1145	853	14	19	56	103	141	166	210	389
ASP*	1040	867	928	15	21	37	49	49	54	59	71
GLN*	41	81	56			57		70	103	194	268
GLU*	1527	1298	1129	38	131	328	475	455	555	820	1770
SER*	475	501	362			187	211	210	275	400	637
GLY*	330	368	416	43	63	166	242	328	352	425	547
ORN*	34	27	26								63
PHE*	47	47	17				7	11	13	16	24
VAL*	140	134	82	25	6	8	12	23	35	49	81
HIS*	26		21								18
LYS*	100	108	35						9	15	19
MET*	31	10								7	20
ILE*	49	40	20					5	8	10	19
LEU*	86	81	31	13	6	6	8	14	21	29	60
THR*	98	86	56			11	20	20	21	26	32
TRP*	11	10	7								2
TYR*	99	88	70	1				5	7	7	12

Amino acid profile: Optimize performance

- In this process, preferred amino acids were the most abundant category; quickly consumed by the 11.5-hour time point to support robust growth
- Intermediately preferred and non-preferred amino acids were lower in abundance, but also consumed early in the ferm
- All amino acids increase as the yeast population growth slows and lyses, while protease continues to produce more amino acids
- Enriching the preferred amino acid category by manipulating nitrogen sources improves performance




Summary: Nitrogen Optimization Strategies

CTE Global provides a comprehensive evaluation of nitrogen in fermentation

- We analyze FAN levels for each contributing source of nitrogen: amino acids, ammonia, and urea
- With our next generation of testing, we will provide amino acid profiling to go beyond FAN to a higher level of performance improvement

Our ultimate goal is to provide our customers with detailed information to **improve their process**



We can help—contact us today.



CTE GLOBAL, INC.

Leading to a brighter tomorrow

cte-global.com

(847) 564-5770