




Development of an ASTM standard for the measurement of “cellulose”

David Mangan
dmangan@neogen.com


Overview

- 1) Megazyme/Neogen and the ASTM process
- 2) The need for a method
- 3) NREL assay and associated biases
- 4) Proposed methods and assay performance
- 5) Next steps
- 6) Conclusion

Megazyme and Polysaccharide Assay Expertise




Development Of An Assay For “Cellulose” Within ASTM



The Need For A Method

The Need For A Method

Example lifecycle of a Renewable Identification Number (RIN)



* RVO = Renewable Volume Obligation

https://www.epa.gov/sites/default/files/2015-08/example_lifecycle_of_a_rin_0.png

The Need For A Method



Fiber Frustration

MAY 10, 2021 BY LISA GIBSON

The stalled U.S. EPA approval of corn kernel fiber-to-ethanol pathways could be having a \$1 billion impact on a 15-billion-gallon-per-year ethanol industry, says Jim Ramm, director of engineering for EcoEngineers. That's assuming 3% of overall production could be from fiber and a \$2 premium. It's worst-case scenario, yes, but it's realistic, nonetheless.



Daily RIN, LCFS & CFP Update

D-Code	Average Price		
	2022	2023	2024
D3	\$3.225	\$3.220	\$2.833
D4	\$0.855	\$0.831	\$0.840
D5	\$0.845	\$0.825	\$0.820
D6	\$0.830	\$0.836	\$0.830

The Need For A Method



Designation: E3181 - 20

Standard Practice for
Determination of the Converted Fraction of Starch and
Cellulosic Content From a Fuel Ethanol Production Facility


$$Ash_{ratio} = (Ash_{AC} / Ash_{BC}) \quad (X5.3)$$

$$CF_c = 1 - ((Cellulosic\ Content_{AC} / Ash_{ratio}) / Cellulosic\ Content_{BC})$$

Where: CF_c = cellulosic converted fraction
 AC = After conversion
 BC = Before conversion

Note:

EPA has defined “cellulosic content” as the sum of **cellulose, hemicellulose and lignin**



NREL Assay and Associated Biases


NREL assay (2021)

Cellulose (2021) 28:1989–2002
<https://doi.org/10.1007/s10570-020-03652-2>


Direct determination of cellulosic glucan content in starch-containing samples

Justin B. Sluiter · Katie P. Michel · Bennett Addison · Yining Zeng · William Michener · Alexander L. Paterson · Frédéric A. Perras · Edward J. Wolfrum

ORIGINAL RESEARCH



Neogen (Modified NREL) Assay (2023)



CHDGG = Cellulose/Hemicellulose-Derived Glucan and Galactan

Bias 1 – Effect Of Yeast In Pre- And Post-fermentation Samples

Cellulose (2021) 28:1989–2002
<https://doi.org/10.1007/s10570-020-03652-2>

ORIGINAL RESEARCH

Direct determination of cellulosic glucan content in starch-containing samples

Justin B. Sluiter · Katie P. Michel · Bennett Addison · Yining Zeng · William Michener · Alexander L. Paterson · Frédéric A. Perras · Edward J. Wolfrum

“We recognize that this method cannot differentiate between β -(1,4) glucans and other β -glucans present. In particular, the measurement of cellulose β -glucans in post-fermentation material will be biased high due to ...fermentation yeast. An attempt to quantify β -(1,3) glucans present in the post fermentation (DDGS) sample using commercially available enzymatic assay (Megazyme K-EBHLG) ... showed significant cross-activity...Until further research can be performed in this area, we are unable to address this bias.”



Designation: E3181 – 20

Standard Practice for
Determination of the Converted Fraction of Starch and
Cellulosic Content From a Fuel Ethanol Production Facility

$$Ash_{ratio} = (Ash_{AC} / Ash_{BC}) \quad (X5.3)$$

$$CF_c = 1 - ((Cellulosic\ Content_{AC} / Ash_{ratio}) / Cellulosic\ Content_{BC})$$


Where:

CF_c = cellulosic converted fraction



AC = After conversion

BC = Before conversion

Bias 1 – Effect Of Yeast In Pre- And Post-fermentation Samples



Bias 1 – Effect Of Yeast In Pre- And Post-fermentation Samples



Sample Description	Treatment	Cellulosic glucan % (w/w) n=2	%CV
Biomass BC	None	1.86	1.59
	+ YDC	1.95	6.80
Biomass AC	None	5.73	4.42
	+ YDC	4.32	2.47

Bias 2 – Loss Of Cellulosic Content Due To NaOH Treatment

Cellulose (2021) 28:1989–2002
<https://doi.org/10.1007/s10570-020-03652-2>

ORIGINAL RESEARCH

Direct determination of cellulosic glucan content in starch-containing samples

Justin B. Sluiter · Katie P. Michel · Bennett Addison · Yining Zeng · William Michener · Alexander L. Paterson · Frédéric A. Perras · Edward J. Wolfrum

“Cellobiose, the β -(1,4)-linked glucan dimer, was chosen as the smallest molecule that retains the linkage indicative of cellulose solubilization; any detectable cellobiose would indicate cellulose loss during enzymatic hydrolysis...we conclude that there is no detectable cellobiose present in the filtrate solutions.”



Carbohydrate Polymers
Volume 51, Issue 3, 15 February 2003, Pages 281-300



Degradation of cellulose under alkaline conditions

Charles J Knill, John F Kennedy



Carbohydrate Polymers
Volume 224, 15 November 2019, 115152



Gelation of cellulose-NaOH solutions in the presence of cellulose fibers

Oona Korhonen^a, Tatiana Budtova^{a b}

Thermal Gelation of Cellulose in a NaOH/Thiourea Aqueous Solution

Lihui Weng, Lina Zhang, Dong Ruan, Lianghe Shi, and Jian Xu

[View Author Information](#)

Cite this: *Langmuir* 2004, 20, 6, 2086–2093
Publication Date: February 18, 2004
<https://doi.org/10.1021/la035995o>
Copyright © 2004 American Chemical Society

Article Views	Altmetric	Citations
1897	3	99

[LEARN ABOUT THESE METRICS](#)



Bioresource Technology
Volume 251, March 2018, Pages 1-6



Cellulase pretreatment for enhancing cold caustic extraction-based separation of hemicelluloses and cellulose from cellulosic fibers

Jianguo Li^{a b c}, Shaokai Zhang^{a c}, Hailong Li^{c d}, Xinhua Ouyang^a, Liulian Huang^a, Yonghao Ni^{a c}, Lihui Chen^a

Review Paper | [Published: 05 November 2015](#)


Cellulose in NaOH–water based solvents: a review

[Tatiana Budtova](#) & [Patrick Navard](#)

Cellulose 23, 5–55 (2016) | [Cite this article](#)

8656 Accesses | 227 Citations | 15 Altmetric | [Metrics](#)

Bias 2 – Loss Of Cellulosic Content Due To NaOH Treatment



Bias 2 – Loss Of Cellulosic Content Due To NaOH Treatment

Sample	n, replicates	% CV	Minus (-) or Plus (+) Ethanolic Precipitation	CHDGG % (w/w)	Recaptured CHDGG % (w/w)
Conventional BC	2	1.5	-	1.69	0.43
	4	2.4	+	2.12	
Conventional AC	2	1.7	-	6.32	0.91
	4	0.3	+	7.23	
CKF Process BC	2	1.3	-	1.82	0.43
	4	0.9	+	2.25	
CKF Process AC	2	0.6	-	4.94	0.79
	4	1.1	+	5.73	

Effect of ethanolic precipitation on the recovery of cellulose in the Neogen assay

Sample	Recovery of CHDGG versus DFO method (%)	
	Without ethanolic precipitation	With ethanolic precipitation
Conventional BC	79.6	95.7
Conventional AC	80.5	102.1
CKF Process BC	88.4	95.8
CKF Process AC	74.4	99.9

Effect of ethanolic precipitation on the recovery of cellulose in the Neogen assay compared to that recovered in the DFO assay

Bias 3 – Exclusion Of Galactose In The Analyte Determined




Figure 3.1 The Leloir pathway. The enzymes of this pathway promote the conversion of β -D-galactose into glucose-1-phosphate which can subsequently be used in glycolysis.



International Review of Cell and Molecular
Biology
Volume 269, 2008, Pages 111-150



Chapter 3 Galactose Metabolism in Yeast— Structure and Regulation of the Leloir Pathway Enzymes and the Genes Encoding Them

Christopher A. Sellick, Robert N. Campbell, Richard J. Reece

BIOTECHNOLOGY
and
BIOENGINEERING

Article

Physiological studies in aerobic batch cultivations of *Saccharomyces cerevisiae* strains harboring the *MEL1* gene

Simon Ostergaard, Christophe Roca, Birgitte Rønnow, Jens Nielsen, Lisbeth Olsson

First published: 31 March 2000 |

[https://doi.org/10.1002/\(SICI\)1097-0290\(20000505\)68:3<252::AID-BIT3>3.0.CO;2-K](https://doi.org/10.1002/(SICI)1097-0290(20000505)68:3<252::AID-BIT3>3.0.CO;2-K) | Citations: 45

Bias 3 – Exclusion Of Galactose In The Analyte Determined

Sample Description	n	%CV	CHDGG		% Cellulosic ethanol	
			% (w/w) DWB		Glucan only	Glucan and galactan
			Glucan only	Glucan and galactan		
Conventional BC	4	2.44	2.21	2.64	-0.01	0.06
Conventional AC	4	0.26	7.63	8.94		
CKF Process BC	4	0.87	2.33	2.79	0.93	1.09
CKF Process AC	4	1.14	6.05	7.25		

Proposed Methods and Performance

Comparison Of Available Open-source Methods

	NREL Method ^a	DFO method ^b	Neogen method ^c	Neogen method ^d
Sample	% Cellulosic ethanol			
Conventional	-0.64	-0.30	0.11	0.08
CKF Process	0.11	0.22	1.21	1.01

^aData from Neogen only

^bCombined data from Neogen, POET, NREL, Novozymes

^cCombined data from Neogen, POET, NREL

^dCombined data from Neogen, POET, NREL - Galactan component excluded

Investigation Into Within And Between Lab Variability

Sample	Lab	n	CHDGG % (w/w)	RSD _r	CHDGG % (w/w)	RSD _R	
Conventional BC	1	13	2.76	2.86	2.77	1.98	
	2	3	2.73	2.15			
	3	3	2.83	10.37			
Conventional AC	1	13	9.47	3.19	9.29	5.76	
	2	3	9.71	0.70			
	3	3	8.69	10.13			
CKF Process BC	1	14	2.87	4.55	2.94	2.25	
	2	3	2.96	3.93			
	3	3	3.00	2.65			
CKF Process AC	1	14	7.71	3.20	7.51	2.87	
	2	3	7.55	2.65			
	3	3	7.28	4.52			
NIST Biomass A	1	14	2.64	8.01	2.50	12.64	
	2	3	2.73	5.33			
	3	3	2.14	4.18			
NIST Biomass B	1	12	6.52	3.88	6.31	2.97	
	2	3	6.21	1.32			
	3	3	6.18	6.67			
				Average RSD_r	4.46	Average RSD_R	4.74

Next Steps

What Could It Mean For The Industry?

On a “per-facility” basis

“Typical” bioethanol production facility ALREADY operating in-situ CKF process

- ✓ Assume **100 Mgal/Yr** capacity
- ✓ Assume **1% cellulosic ethanol** based on Neogen methodology
- ✓ Assume **\$1.60** D6-D3 RIN spread

Predicts ~\$1.6m annual benefit with no operational/manufacturing modification before:

- ✓ YDC license fee
- ✓ Analytic lab fees
- ✓ Additional audit costs

What Could It Mean For The Industry?

Using LCFS to approximate the quantity of CKF processes in situ today

Phase 1 – Facilities that have already adopted cellulase-assisted fermentation technology:

- ✓ Ethanol production capacity approved in the LCFS program is **~5.8Bgal**
- ✓ Assume **1% cellulosic ethanol** based on Neogen methodology
- ✓ Assume **\$1.60** D6-D3 RIN spread


Predicts ~\$93m uplift for those plants on a no-change basis

Phase 2 – Facilities that currently utilize conventional fermentation technology:

- ✓ Ethanol production capacity NOT approved in the LCFS program is **~10.5Bgal**
- ✓ Same assumptions as above

Predicts potential ~\$168m D3 RIN “carrot”

How Do We Get There?



How Will It Work?

Phase 1

Producer objective:
Prepare for facility
registration by EPA

E-RDYDC



Tasks:

- Purchase E-RDYDC under supply agreement and perform in-house analysis or outsource to your analytical partner
- Gather % cellulosic ethanol data from multiple fermentations
- Plant fermentation process development (if required)
- Predict % cellulosic ethanol

<https://www.epa.gov/fuels-registration-reporting-and-compliance-help/how-register-new-renewable-fuel-producer-renewable>

Phase 2

Producer objective:
Generate D3 RINs

E-YDC



Tasks:

- Order YDC (\$FOC) under supply agreement and perform in-house analysis or outsource to your analytical partner
- Generate % cellulosic ethanol data from multiple fermentations
- Complete 3rd party engineering audit and submit to EPA for facility registration
- Re-submit data to EPA every 500K gallons for re-certification
- Pay a license fee to Neogen based on D3 \$ premium generated on a monthly basis

<https://www.govinfo.gov/content/pkg/CFR-2022-title40-vol19/pdf/CFR-2022-title40-vol19-sec80-1451.pdf> - see page 4

Summary

- A **VCSB method** has been developed for the measurement of the relevant portion of **cellulosic content** as defined by the EPA
- This method is currently under consideration for adoption as a standard within **ASTM**
- Once approval has been finalized, **EPA approval** is expected to quickly follow based on existing guidance
- For those plants already generating cellulosic ethanol, conversion of the relevant D6 to D3 RINs will be possible with no operational changes required, and will deliver significant financial benefit
- For those plants not currently generating cellulosic ethanol, could this be the additional benefit needed to justify the change?!

For more information on how to join ASTM and help support the VCSB method approval, contact David Mangan (dmangan@neogen.com)

For more information on the commercial aspects of accessing the analytical technology, contact Matt Nichols (mnichols@neogen.com)

Acknowledgements

Neogen

Ruth Ivory
Anna Draga
Tadas Kargelis
Amaya O’Cochláin
Lucie Charmier
Artur Rogowski
Andrea Mascherpa

NREL

Justin Sluiter
Katie Michel

POET

Bart Plocher
Melissa Tille

Novozymes

David Gogerty
Geoff Moxley

TCD

John O’Brien
Manuel Ruether

CCRC

Paristoo Azadi

ASTM standard collaborators

Too many to list!

Thank you
for listening



Any questions?