



Where science
& creativity meet

INTRODUCTION TO IFF

Keaton Albers– Senior Applications Associate II

IFF - APPLIED INNOVATION CENTER

Cedar Rapids, Iowa

Lab-based plant support

- Pre-trial testing
- Trial evaluation
- Optimization
- Troubleshooting



Liquefaction Services

- Cook studies
- Solubility
- Cations (Sodium, etc.)



Fermentation Services

- Prop and ferm studies
- DP4+ composition
- Detailed sugar analysis
- HPLC checks
- Residual starch
- Nitrogen measurements
- Inhibitors (fusels, sodium, sulfite, organic acids, etc.)





XCELIS® Ethanol Solutions

HPLC TROUBLESHOOTING AND TRAINING

HPLC TROUBLESHOOTING

Diagnostic evaluation

Two main types of problems



Hardware issues

- Column or other consumables
- Major system error



Software issues

- Integration or calibration issues



How is this different from normal operations?

- Keep records of normal operations (or be able to find previous data)
- Run check standards

Start with hardware problems

- Quick to check

HARDWARE PROBLEMS



Changes in physical parameters usually indicate a hardware issue



Quick and easy-to-spot changes

- Initial check can be done in <5 minutes with good record keeping



What changes/observations are the operators/lab manager seeing when the HPLC is operating?



How is this different from normal operations?

Keep records of normal operations (or be able to find previous data)

- Pressure
- Retention time
- Peak shape (especially width and tailing)
- Area (choose a Standard)

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 1:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 1:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 1:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 1:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 1:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 2:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS




High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 2:

Observed pressure:



Replace plugged part

HARDWARE PROBLEMS



High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 2:



Replace plugged part

HARDWARE PROBLEMS



High pressure



First step: Find the pressure source

- Start from detector, work backward, and remove components until pressure source is found
- Example 2:



Replace plugged part

HARDWARE PROBLEMS

Old Column

Retention time:


Compare retention times to previous runs

- Standards work best
- If retention times changed, check peak windows and peak IDs

Peak shape

- Wide peaks and tailing peaks both usually indicate column issues
- Often only solution is to replace column
- Can sometimes indicate old tubing or pump issues

Blue column needs replaced



HARDWARE PROBLEMS


Changing Columns

Retention Time Changes

Peak IDs can change as small retention time changes can move peaks outside of peak ID window

Coelution changes

Peak splitting, especially around DP4+, DP2, DP1



HARDWARE PROBLEMS

Changing Columns

Retention Time Changes


Peak IDs can change as small retention time changes can move peaks outside of peak ID window

Coelution changes

Peak splitting, especially around DP4+, DP2, DP1

Improving column performance can cause problems

Same issues as poor performance, but in reverse



HARDWARE PROBLEMS

Inconsistent Areas



Standard area check



Compare area of standard or check standards to previous runs of the same standard



Change in area indicates a possible HPLC problem




If no other issues are present (integration, peak shape, etc.), changes in area may indicate a problem with the standard, autosampler, or detector

| Date | Pressure | DP1 RT | DP1 Area | EtOH RT | EtOH Area |
|------|----------|--------|----------|---------|-----------|
| 4/3 | 466 | 9.82 | 100563 | 20.1 | 100468 |
| 4/5 | 469 | 9.83 | 101571 | 19.8 | 99731 |
| 4/10 | 482 | 9.81 | 103648 | 19.9 | 99464 |
| 4/12 | 470 | 9.85 | 99615 | 19.9 | 100753 |
| 4/14 | 476 | 9.83 | 80640 | 20.0 | 82376 |

REPRODUCIBILITY CHECK

- ✓ Make Standards and check standards and place in 5-10 vials



- ✓ Inject the vials once or twice per day as follows (vary time of day if once/day) on each HPLC :

| | | | | | | |
|----------|----------------|----------------|-----------------|-------------------|-------------------|-------------------|
| Day 1 AM | Low Standard 1 | Med Standard 1 | High Standard 1 | Check Standard 1 | Check Standard 1 | Check Standard 1 |
| Day 1 PM | | | | Check Standard 2 | Check Standard 2 | Check Standard 2 |
| Day 2 AM | Low Standard 2 | Med Standard 2 | High Standard 2 | Check Standard 3 | Check Standard 3 | Check Standard 3 |
| Day 2 PM | | | | Check Standard 4 | Check Standard 4 | Check Standard 4 |
| Day 3 AM | Low Standard 3 | Med Standard 3 | High Standard 3 | Check Standard 5 | Check Standard 5 | Check Standard 5 |
| Day 3 PM | | | | Check Standard 6 | Check Standard 6 | Check Standard 6 |
| Day 4 AM | Low Standard 4 | Med Standard 4 | High Standard 4 | Check Standard 7 | Check Standard 7 | Check Standard 7 |
| Day 4 PM | | | | Check Standard 8 | Check Standard 8 | Check Standard 8 |
| Day 5 AM | Low Standard 5 | Med Standard 5 | High Standard 5 | Check Standard 9 | Check Standard 9 | Check Standard 9 |
| Day 5 PM | | | | Check Standard 10 | Check Standard 10 | Check Standard 10 |

- ✓ Results (%wt/v) should be the same for each check standard run
- ✓ If results vary, HPLC system has a reproducibility problem
- ✓ This can prove that your system has an issue

SOFTWARE PROBLEMS




If there is no obvious instrument malfunctions, check data analysis problems




Common data analysis problems

- Peak identifications
- Peak integrations
- Calibration errors

STANDARD CORN-TO-ETHANOL FERMENTATION HPLC CHROMATOGRAM




STANDARD CORN-TO-ETHANOL FERMENTATION HPLC CHROMATOGRAM




DP4+ INTEGRATION

- ✓ DP4+ peak is sometimes split
- ✓ Peak splitting can be variable
- ✓ Peak splitting can change
 - Especially with phytase treatment
 - Column changes
- ✓ Best practice is to be consistent
- ✓ Peak groups are helpful in integrating together



DP2 INTEGRATION

- ✓ DP2 and phosphate coelute
- ✓ Changes in DP2 or phosphate can change the way the peaks are integrated
- ✓ Peak IDs can shift with changing peak sizes



DP2 INTEGRATION OPTIONS



For proper splitting

- Shoulder sensitivity
- Overall sensitivity (threshold/smoothing)
- Apply over set time window if possible




For proper peak IDs



- Change peak windows
- Add phosphate peak



Changes may adversely affect peaks at other time points





OTHER INTEGRATION DIFFICULTIES- BASELINE CHANGES




| Baseline type: | Valley-Valley | Flat |
|----------------|---------------|-------|
| Batch: | 4672 | 4673 |
| DP4+ | 0.42 | 0.49 |
| DP3 | 0.06 | 0.01 |
| DP2 | 0.26 | 0.36 |
| DP1 | 0.06 | 0.14 |
| Total sugars | 0.8 | 1 |
| Lactic | 0.06 | 0.07 |
| Glycerol | 0.75 | 0.72 |
| Acetic Acid | 0.08 | 0.13 |
| Ethanol | 14.63 | 14.51 |

OTHER INTEGRATION DIFFICULTIES- PEAK SHOULDERS



FUSEL INTERFERENCE IN ACETATE PEAK

Plant was experiencing sluggish fermentations and seeing unusually high acetate in ferm samples. AIC fusel analysis confirmed fusel recycle issue.



Fusels from previous injections carried over into next run




Can elute near acetate and interfere with acetate quantitation

FUSEL INTERFERENCE IN ACETATE PEAK

✓ Black: Liquefact injected first run of the day

✓ Blue: Liquefact injected after drop sample with high fusels



BAD CALIBRATION CURVE

✓ Check Standard

✓ Compare areas

- Changes in areas should roughly match changes in amounts


✓ If areas are similar but amounts are different, there may be a calibration curve issue

✓ In this case, DP4+ was not integrated properly in the Standard

| | Batch | Area | Amount | Area %Δ | Amount %Δ |
|----------|-------|---------|--------|---------|-----------|
| DP4+ | 5000 | 842242 | 0.433 | 0 | 0 |
| | 5050 | 1195681 | 0.607 | 42 | 40 |
| | 5075 | 1064758 | 0.797 | 26 | 84 |
| | 5100 | 1224226 | 0.933 | 45 | 115 |
| DP1 | 5000 | 112252 | 0.084 | 0 | 0 |
| | 5050 | 53686 | 0.039 | -52 | -54 |
| | 5075 | 50499 | 0.038 | -55 | -55 |
| | 5100 | 143093 | 0.114 | 27 | 36 |
| Glycerol | 5000 | 1746929 | 1.569 | 0 | 0 |
| | 5050 | 1685821 | 1.515 | -3 | -3 |
| | 5075 | 1811692 | 1.629 | 4 | 4 |
| | 5100 | 1809201 | 1.625 | 4 | 4 |


BAD CALIBRATION CURVE- ETHANOL

- ✓ Unexpected or drifting results
- ✓ Sudden change in EtOH values
 - Especially after a standard change
- ✓ Old standards lose ethanol, but software's calibration level stays the same
- ✓ 5% decrease in high standard area results in sample [EtOH] changing from 15.06% to 15.78%




BAD CALIBRATION CURVE- ETHANOL

- ✓ Unexpected or drifting results
- ✓ Sudden change in EtOH values
 - Especially after a standard change
- ✓ Old standards lose ethanol, but software's calibration level stays the same
- ✓ 5% decrease in high standard area results in sample [EtOH] changing from 15.06% to 15.78%




BAD CALIBRATION CURVE- ETHANOL

- ✓ Unexpected or drifting results
- ✓ Sudden change in EtOH values
 - Especially after a standard change
- ✓ Old standards lose ethanol, but software's calibration level stays the same
- ✓ 5% decrease in high standard area results in sample [EtOH] changing from 15.06% to 15.78%




BAD CALIBRATION CURVE- ETHANOL

- ✓ Unexpected or drifting results
- ✓ Sudden change in EtOH values
 - Especially after a standard change
- ✓ Old standards lose ethanol, but software's calibration level stays the same
- ✓ 5% decrease in high standard area results in sample [EtOH] changing from 15.06% to 15.78%



BAD CALIBRATION CURVE- ETHANOL


- ✓ Unexpected or drifting results
- ✓ Sudden change in EtOH values
 - Especially after a standard change
- ✓ Old standards lose ethanol, but software's calibration level stays the same
- ✓ 5% decrease in high standard area results in sample [EtOH] changing from 15.06% to 15.78%



BAD CALIBRATION CURVE- ETHANOL

Resulting Data

✓ Changes in amounts do not correlate with changes in peak area



HPLC PROBLEM CHECKLIST



Identify what is different



Before an issue arises: know what is normal

- Run frequent check standards



After the issue

- Check for physical problems
- Check for integration/peak ID issues
- Check for calibration issues
- Check for reproducibility



If still having issues- note timelines and what could have changed



THANK YOU!

Questions?

STAY CONNECTED



www.xcelis.com



<https://www.linkedin.com/showcase/xcelis-ethanol-solutions/>

Grain Changers Community



Join Grain Changers for exclusive content for our valued partners.

Register for free at:

www.xcelis.com/grain-changers/

©2023 International Flavors & Fragrances Inc. (IFF). IFF, the IFF Logo, and all trademarks and service marks denoted with ™, SM or ® are owned by IFF or affiliates of IFF unless otherwise noted. The information provided herein is based on data IFF believes, to the best of its knowledge, reliable and applies only to the specific material designated herein as sold by IFF. The information contained herein does not apply to use of the material designated herein in any process or in combination with any other material and is provided at the request of and without charge to our customers. Accordingly, IFF cannot guarantee or warrant such information and assumes no liability for its use. Other than as may be expressly set forth in a contract of sale, IFF makes no warranty, express or implied, as to the material set forth herein, including the warranty of merchantability or fitness for a particular use.

