Save on PYREX® Glassware

Exclusive Promotional Offer:

Purchase 3 cases of PYREX or PYREXPLUS® Reusable Glassware and Get 1 Case of Selected PYREX Glassware at No Additional Cost

For purchases April 15, 2024 through June 15, 2024



PYREX® Beakers, Standard, Low Form, Griffin

Corning Cat. No.	Description	Qty/Pk	Qty/Cs
1000-50	50 mL beaker, graduated	12	48
1000-100	100 mL beaker, graduated	12	48
1000-150	150 mL beaker, graduated	12	48



PYREX® Wide Mouth Erlenmeyer Flasks, with GL45 Screw Cap

Cat. No.	Description	Qty/Pk	Qty/Pk Qty/Cs	
4995-250	250ML Wide Mouth Erlenmeyer Flask	6	6	
4995-500	500ML Wide Mouth Erlenmeyer Flask	6	6	
4995-1L	1L Wide Mouth Erlenmeyer Flask	6	6	
4995-2L	2L Wide Mouth Erlenmeyer Flask	1	1	



PYREX[®] Cylinders, Single Metric Scale, Graduated Cylinder, TC

Corning Cat. No.	Description	Qty/Pk	Qty/Pk Qty/Cs	
3022-10	10 mL cylinder, with funnel top	1	24	
3022-50	50 mL cylinder	1	18	
3022-100	100 mL cylinder	1	12	
3022-250	250 mL cylinder	1	12	
3022-500	500 mL cylinder	1	8	
3022-1L	1L cylinder	1	1	



PYREX[®] Volumetric Flask, Class A, Heavy Duty, with Glass Standard Taper Stopper

Cat. No.	Description		Qty/Pk Qty/Cs	
5635-10	10 ml Volumetric Flask	6	6	
5635-50	50 ml Volumetric Flask	6	6	
5635-100	100 ml Volumetric Flask	6	6	
5640-25	25 ml Volumetric Flask	6	12	
5640-50	50 ml Volumetric Flask	6	12	
5640-100	100 ml Volumetric Flask	6	12	
5640-250	250ml Volumetric Flask	6	12	
5640-500	500ml Volumetric Flask	6	12	
5640-1L	1L Volumetric Flask	1	6	
5642-100	100 ml Volumetric Flask	6	12	
5840-200	200 mL Phosphoric Acid Vol Flasks	6	12	
5840-250	250 mL Phosphoric Acid Vol Flasks	6	12	



PYREX® Heavy Wall Filtering Flasks with Sidearm Tubulation

Corning Cat. No.	Description	Qty/Pk Qty/Cs	
5360-125	125 mL Micro Filtering Flasks with Sidearm Tubulation	6	18
5340-250	250 mL Heavy Wall Filtering Flasks with Sidearm Tubulation	6	18
5340-1L	1L Heavy Wall Filtering Flasks with Sidearm Tubulation	6	12
5340-2L	2L Heavy Wall Filtering Flasks with Sidearm Tubulation	1	12
5340-4L	4L Heavy Wall Filtering Flasks with Sidearm Tubulation	1	6



To redeem this promotion, please reach out to your local Corning Account Manager:

laboratory environments, making it the indispensable partner for chemists and research scientists for more than 100 years.

Name



Exceptional strength and stability. Brilliant transparency. Superior resistance to chemicals, contaminants, and drastic changes in temperature. These are among the key properties that enable PYREX® glass to endure even the toughest

Email

Phone

Terms and Conditions: Offer valid in the U.S. and Canada only and void where prohibited by government, law, or company policy. Offer valid for end-user customers only for purchases April 15, 2024 through June 15, 2024 or while supply lasts. Redemption information must be received by July 1, 2024. Please allow 4-6 weeks for delivery. For additional product or technical information, visit www.corning.com/life-sciences, call 1.800.492.1110, or contact your Corning account manager. Free product must be of equal or lesser value of purchased product.